

March 2, 2026

Laboratory 02

Invite laboratory members (PI only)

Alternative way

19th IHIWS database system development team

Secret registration token (PI only)

The screenshot shows a web browser window with the URL <https://data.ihiw19.org>. The page title is "IHIW19 Dashboard" and the user is logged in as "IHIW19 PIDEMO1". The main heading is "Dashboard" with a sub-heading "Welcome to IHIW19 Management System, IHIW19 PIDEMO1".

The "Metadata Completion Status" section contains two panels:

- HLA Metadata** (blue header):
 - Sample Metadata (Stage 1): Sample: 0 of 0 completed (0%)
 - Subject Metadata (Stage 2): Subject: 0 of 0 completed (0%)
 - No HLA metadata found
- Antibody Metadata** (yellow header):
 - Sample Metadata (Stage 1): Sample: 0 of 0 completed (0%)
 - Subject Metadata (Stage 2): Subject: 0 of 0 completed (0%)
 - No Antibody metadata found

At the bottom, there are two navigation cards:

- Lab View**: Manage lab members, upload files, and configure sharing settings with projects. Includes a "Go to Lab View" button.
- Project View**: Manage projects, view shared files, and import/export data. Includes a "Go to Project View" button.









A callout box with a green border points to the "Go to Lab View" button, containing the text: "Click 'Go to Lab View'".

© 2025 IHIW19 Management System

Add Lab Member (PI only)

The screenshot shows a web browser window with the URL <https://data.ihw19.org/lab>. The page title is "Lab View" and the lab name is "PRDISH". The user is logged in as "Manabu Production PI Ishii".

On the left sidebar, the "Lab Members" menu item is highlighted. The main content area shows a "Lab Members" section with a table of members and an "Add Lab Member" button. A callout box points to the button with the text "Click 'Add Lab Member'".

Name	Email	Role	Status	Action
Manabu Production PI Ishii	[blurred]	PI	Active	
Manabu Production LAB Member Ishii	[blurred]	Lab Member	Active	 
Manabu LabMember 2 Ishii	[blurred]	PI	Active	 
Manabu Production ADMIN Ishii	[blurred]	PI	Active	 
Manabu LabMember3	[blurred]	PI	Active	 

Add Lab Member (PI only)

The screenshot shows the IHIW19 Lab View interface. A modal dialog titled "Add Lab Member" is open, containing three input fields: "First Name *", "Last Name *", and "Email *". Below the fields are "Cancel" and "Create" buttons. Three callout boxes on the left point to the input fields with the text "Fill 'First Name'", "Fill 'Last Name'", and "Fill 'Email'". A fourth callout box on the right points to the "Create" button with the text "Click 'Create' Button".

Lab View | IHIW19 Management

https://data.ihw19.org/lab

IHIW19 User Guide Dashboard Lab View Project View

Manabu Production PI Ishii

Add Lab Member

First Name *

Last Name *

Email *

Cancel Create

Fill "First Name"

Fill "Last Name"

Fill "Email"

Click "Create" Button

Name	Role	Status	Actions
Manabu LabMember 1 Ishii	Lab Member	Active	[Edit] [Key]
Manabu LabMember 2 Ishii	PI	Active	[Edit] [Key]
Manabu Production ADMIN Ishii	PI	Active	[Edit] [Key]
Manabu LabMember3	PI	Active	[Edit] [Key]

Add Lab Member (PI only)

Lab View | IHIW19 Management

https://data.ihw19.org/lab

IHIW19 User Guide Dashboard Lab View Project View

Manabu Production PI Ishii

Lab View

Lab Members

- Project Linkage
- Genotype Data
- Antibody Testing Data
- Sample Metadata
- Subject Metadata
- Token Management

Lab Member Created

Please share the following password securely with the new lab member. They will be required to change it on first login.

Generated Password:

[Redacted Password]

Close

Copy "Generated Password"

Push "Close" Button

Name	Role	Status	Action
Manabu Production LAB Member Ishii	Lab Member	Active	[Add] [Key]
Manabu LabMember 2 Ishii	PI	Active	[Add] [Key]
Manabu LabMember3 Ishii	PI	Active	[Add] [Key]
Manabu LabMember4 Ishii	Lab Member	Active	[Add] [Key]

Provide the email address and the generated password (Initial Password) to the Lab Member.