#### User Guide for the 19th IHIWS Database System

September 23, 2025

## Laboratory 07

# Upload Antibody Testing Data

One Lambda


Werfen


19th IHIWS database system development team


#### **Antibody Testing Data Format**


Our database system accepts antibody testing data in Comma-Separated Values (CSV) format exported from **One Lambda** and **Werfen (Immucor)** software.


User guides for exporting CSV files from each software application are currently in preparation.














https://data.ihiw19.org/lab/#antibody-testing-data 🙎 ゲスト IHIW19 Tsuyoshi Hachiya2 **Import Antibody Data** X CSV File head-IWSampleCl2025-05-13.csv ファイルを選択 Lab: GAJ Upload a CSV file in OneLambda or Immucor format. Lab Members **m** Delete Selected **Typing Method Information** Project Linkag Test Type Vendor ▼ Genotype Data One Lambda ThermoFisher SAB Select "Test Type" Select "Vendor" Software Instrument ■ Sample Metada Select Select "Software" Luminex FM3D **Fusion** "Instrument" Software version Additional Information Token Manage Input "Software Input "Additional e.g. 1.0 Optional version" Information" **Upload** Close Click "Upload"







ILIW10 Dachboard Lab View

Individual Antibody Details: GTI 64

X

Sample ID: GTI 64

Lab Code: JPNHAC2

Individual Specificity/MFI Data:

Specificity	Raw MFI	Normalized MFI	Bead Classification	Bead Rank
A*01:01	243	5	Negative	2
A*02:01	141	3	Negative	2
A*02:02	158	3	Negative	2
A*02:03	194	4	Negative	2
A*02:05	155	3	Negative	2
A*03:01	158	3	Negative	2
A*11:01	139	3	Negative	2
A*11:02	165	3	Negative	2
A*23:01	782	15	Positive	8

Total records: 9 | Positive: 1 | Negative: 8

Data detail will be shown

Close

