

March 4, 2026

Laboratory 10

Data Visibility

Excel

19th IHIWS database system development team

Manage Data Visibility

https://data.ihw19.org

IHW19 [User Guide](#) **Dashboard** [Lab View](#) [Project View](#)

 Test User ▾

Dashboard

Welcome to IHW19 Management System, Test User

Metadata Completion Status

HLA Metadata


Sample Metadata (Stage 1)

Sample: 0 of 0 completed (0%)



Subject Metadata (Stage 2)

Subject: 0 of 0 completed (0%)




No HLA metadata found

Antibody Metadata

Sample Metadata (Stage 1)

Sample: 0 of 0 completed (0%)



Subject Metadata (Stage 2)

Subject: 0 of 0 completed (0%)



No Antibody metadata found

Click "Lab View"

Lab View

Manage lab members, upload files, and configure sharing settings with projects.

[Go to Lab View](#)

Project View

Manage projects, view shared files, and import/export data.

[Go to Project View](#)

Manage Data Visibility

https://data.ihw19.org/lab

IHW19 User Guide Dashboard **Lab View** Project View

Test User ▾

Lab View

Lab: DOCLAB01

- Lab Members**
- Project Linkage
- Genotype Data
- Antibody Testing Data
- Sample Metadata
- Subject Metadata
- Data Visibility**
- Token Management

Lab Members

[+ Add Lab Member](#)

Name	Email	Role	Status	Action
Test User	testuser@example.com	PI	Active	

Click "Data Visibility"

Manage Data Visibility

https://data.ihw19.org/lab

IHW19 User Guide Dashboard **Lab View** Project View

Test User ▾

Lab View

Lab: DOCLAB01

- Lab Members
- Project Linkage
- Genotype Data
- Antibody Testing Data
- Sample Metadata
- Subject Metadata
- Data Visibility**
- Token Management

Data Visibility

[Download Excel](#)

[Upload Excel](#)

Control which subjects' data is visible to each linked project. Changes take effect immediately for project leaders viewing ProjectView.

Subject_ID	System_Subject_ID	HLA Typing Workshop 2025	Antibody Screening 2025
SUBJ_001	IHW19SUB0000000001	YES	YES
SUBJ_002	IHW19SUB0000000002	YES	YES
SUBJ_003	IHW19SUB0000000003	YES	YES
SUBJ_004	IHW19SUB0000000004	YES	YES

Data Visibility matrix shows per-subject, per-project visibility settings

Manage Data Visibility

https://data.ihw19.org/lab

IHW19 User Guide Dashboard **Lab View** Project View

Test User ▾

Lab View

Lab: DOCLAB01

- Lab Members
- Project Linkage
- Genotype Data
- Antibody Testing Data
- Sample Metadata
- Subject Metadata
- Data Visibility**
- Token Management

Data Visibility

[Download Excel](#)

[Upload Excel](#)

Control which subjects' data is visible to each linked project. Changes take effect immediately for project leaders viewing ProjectView.

Subject_ID	System_Subject_ID	HLA Typing Workshop 2025	Antibody Screening 2025
SUBJ_001	IHW19SUB0000000001	YES	YES
SUBJ_002	IHW19SUB0000000002	YES	YES
SUBJ_003	IHW19SUB0000000003	YES	YES
SUBJ_004	IHW19SUB0000000004	YES	YES

Subject_ID: Original subject identifier from your lab

YES (green) = visible to project
NO (red) = hidden from project

Manage Data Visibility

https://data.ihw19.org/lab

Lab View

Lab: DOCLAB01

- Lab Members
- Project Linkage
- Genotype Data
- Antibody Testing Data
- Sample Metadata
- Subject Metadata
- Data Visibility**
- Token Management

Click "Download Excel"

[Download Excel](#) [Upload Excel](#)

Data Visibility

Control which subjects' data is visible to each linked project. Changes take effect immediately for project leaders viewing ProjectView.

Subject_ID	System_Subject_ID	HLA Typing Workshop 2025	Antibody Screening 2025
SUBJ_001	IHW19SUB0000000001	YES	YES
SUBJ_002	IHW19SUB0000000002	YES	YES
SUBJ_003	IHW19SUB0000000003	YES	YES
SUBJ_004	IHW19SUB0000000004	YES	YES

Manage Data Visibility

https://data.ihw19.org/lab

IHIW19 User Guide Dashboard **Lab View** Project View

Test User ▾

Lab View

Lab: DOCLAB01

- Lab Members
- Project Linkage
- Genotype Data
- Antibody Testing Data
- Sample Metadata
- Subject Metadata
- Data Visibility**
- Token Management

Excel file downloaded successfully

Data Visibility

Download Excel

Upload Excel

Visibility Excel downloaded successfully!

Control which subjects' data is visible to each linked project. Changes take effect immediately for project leaders viewing ProjectView.

Subject_ID	System_Subject_ID	HLA Typing Workshop 2025	Antibody Screening 2025
SUBJ_001	IHIW19SUB00000000001	YES	YES
SUBJ_002	IHIW19SUB00000000002	YES	YES
SUBJ_003	IHIW19SUB00000000003	YES	YES
SUBJ_004	IHIW19SUB00000000004	YES	YES

Manage Data Visibility

Open the downloaded Excel file

	A	B	C	D	E
1	Subject_ID	System_Subject_ID	HLA Typing Workshop	Antibody Screening 2025	
2	SUBJ_001	IHIW19SUB0000000001	YES	YES	
3	SUBJ_002	IHIW19SUB0000000002	YES	YES	
4	SUBJ_003	IHIW19SUB0000000003	YES	YES	
5	SUBJ_004	IHIW19SUB0000000004	YES	YES	

“Subject_ID” and “System_Subject_ID” columns are read-only (gray background)

Manage Data Visibility

	A	B	C	D	E
1	Subject_ID	System_Subject_ID	HLA Typing Workshop	Antibody Screening 2025	
2	SUBJ_001	IHIW19SUB0000000001	NO	NO	
3	SUBJ_002	IHIW19SUB0000000002	YES	YES	
4	SUBJ_003	IHIW19SUB0000000003	YES	NO	
5	SUBJ_004	IHIW19SUB0000000004	YES	YES	

Project columns: Select YES or NO
YES = data visible to project leader
NO = data hidden from project leader

Manage Data Visibility

	A	B	C	D	E
1	Subject_ID	System_Subject_ID	HLA Typing Workshop	Antibody Screening 2025	
2	SUBJ_001	IHIW19SUB0000000001	NO	NO	
3	SUBJ_002	IHIW19SUB0000000002	YES	YES	
4	SUBJ_003	IHIW19SUB0000000003	YES	NO	
5	SUBJ_004	IHIW19SUB0000000004	YES	YES	



Change YES to NO (or NO to YES) for each subject and project

Only YES and NO values are accepted. Other values will cause an error during upload.

Manage Data Visibility

	A	B	C	D	E
1	Subject_ID	System_Subject_ID	HLA Typing Workshop	Antibody Screening 2025	
2	SUBJ_001	IHIW19SUB0000000001	NO	NO	
3	SUBJ_002	IHIW19SUB0000000002	YES	YES	
4	SUBJ_003	IHIW19SUB0000000003	YES	NO	
5	SUBJ_004	IHIW19SUB0000000004	YES	YES	



Save the updated Excel file

Manage Data Visibility

https://data.ihw19.org/lab

IHIW19 User Guide Dashboard **Lab View** Project View

Test User ▾

Lab View

Lab: DOCLAB01

Click "Upload Excel"

Download Excel

Upload Excel

Data Visibility

Visibility Excel downloaded successfully!

Control which subjects' data is visible to each linked project. Changes take effect immediately for project leaders viewing ProjectView.

Subject_ID	System_Subject_ID	HLA Typing Workshop 2025	Antibody Screening 2025
SUBJ_001	IHIW19SUB00000000001	YES	YES
SUBJ_002	IHIW19SUB00000000002	YES	YES
SUBJ_003	IHIW19SUB00000000003	YES	YES
SUBJ_004	IHIW19SUB00000000004	YES	YES

Select the updated Excel file and click "Open"

- Lab Members
- Project Linkage
- Genotype Data
- Antibody Testing Data
- Sample Metadata
- Subject Metadata
- Data Visibility**
- Token Management

Manage Data Visibility

https://data.ihw19.org/lab

IHW19 User Guide Dashboard **Lab View** Project View

Test User ▾

Lab View

Lab: DOCLAB01

- Lab Members
- Project Linkage
- Genotype Data
- Antibody Testing Data
- Sample Metadata
- Subject Metadata
- Data Visibility**
- Token Management

Data Visibility

Download Excel

Upload Excel

Visibility Excel uploaded successfully!

Control which subjects' data is visible to each linked project. Changes take effect immediately for project leaders viewing ProjectView.

Subject_ID	System_Subject_ID	HLA Typing Workshop 2025	Antibody Screening 2025
SUBJ_001	IHIW19SUB00000000001	NO	NO
SUBJ_002	IHIW19SUB00000000002	YES	YES
SUBJ_003	IHIW19SUB00000000003	YES	NO
SUBJ_004	IHIW19SUB00000000004	YES	YES

Upload completed successfully. The visibility matrix is updated.

Manage Data Visibility

https://data.ihw19.org/lab

IHIW19 User Guide Dashboard **Lab View** Project View

Test User ▾

Lab View

Lab: DOCLAB01

👤 Lab Members
🔗 Project Linkage
🕒 Genotype Data
🔪 Antibody Testing Data
🏷️ Sample Metadata
👤 Subject Metadata
👁️ Data Visibility
🔑 Token Management

Data Visibility

Download Excel

Upload Excel

Visibility Excel uploaded successfully!

Control which subjects' data is visible to each linked project. Changes take effect immediately for project leaders viewing ProjectView.

Subject_ID	System_Subject_ID	HLA Typing Workshop 2025	Antibody Screening 2025
SUBJ_001	IHIW19SUB00000000001	NO	NO
SUBJ_002	IHIW19SUB00000000002	YES	YES
SUBJ_003	IHIW19SUB00000000003	YES	NO
SUBJ_004	IHIW19SUB00000000004	YES	YES

Visibility settings are now updated for the selected subjects

Changes take effect immediately. Project leaders will see or not see the data based on the updated settings.