

May 22, 2026

Laboratory 11

Genotype Sample Visibility

Excel

19th IHIWS database system development team

Manage Genotype Sample Visibility

<https://data.ihw19.org>

IHW19

[User Guide](#)

Dashboard

[Lab View](#)

[Project View](#)

 Test User ▾

Dashboard

Click "Lab View"

Welcome to IHW19 Management System, Test User

Metadata Completion Status

HLA Metadata

Sample Metadata (Stage 1)

Sample: 0 of 4 completed (0%)

Subject Metadata (Stage 2)

Subject: 0 of 4 completed (0%)

No HLA metadata found

Antibody Metadata

Sample Metadata (Stage 1)

Sample: 0 of 0 completed (0%)

Subject Metadata (Stage 2)

Subject: 0 of 0 completed (0%)

No Antibody metadata found

Lab View

Manage lab members, upload files, and configure sharing settings with projects.

[Go to Lab View](#)

Project View

Manage projects, view shared files, and import/export data.

[Go to Project View](#)

Manage Genotype Sample Visibility

https://data.ihw19.org/lab

IHW19 User Guide Dashboard **Lab View** Project View

Test User ▾

Lab View

Lab: DOCLAB02

Samples not visible to any project: HLA: 0 / Antibody: 0

Samples visible per project:

- HLA Typing Workshop 2025: HLA 4 (hidden: 0) / Antibody 0 (hidden: 0)
- Antibody Screening 2025: HLA 4 (hidden: 0) / Antibody 0 (hidden: 0)

Lab Members

Project Linkage

Genotype Data

Antibody Testing Data

Sample Metadata

Subject Metadata



Data Visibility

Genotype Sample Visibility

Token Management

Lab Members

Add Lab Member

Name	Email	Role	Status	Action
Test User	testuser@example.com	PI	Active	
	testuser_metadata@example.com	Lab Member	Active	 

Click "Genotype Sample Visibility"

Manage Genotype Sample Visibility

<https://data.ihw19.org/lab/#genotype-sample-visibility>

IHIW19 [User Guide](#) [Dashboard](#) **[Lab View](#)** [Project View](#)

 Test User ▾

Lab View

Lab: DOCLAB02


Samples not visible to any project: HLA: 0 / Antibody: 0

Samples visible per project:

- HLA Typing Workshop 2025: HLA 4 (hidden: 0) / Antibody 0 (hidden: 0)
- Antibody Screening 2025: HLA 4 (hidden: 0) / Antibody 0 (hidden: 0)

 Lab Members

 Project Linkage

 Genotype Data

 Antibody Testing Data

 Sample Metadata

 Subject Metadata

 Data Visibility

 **Genotype Sample Visibility**

 Token Management

Genotype Sample Visibility

 Download Excel

 Upload Excel

Control which individual Genotype samples are visible to each linked project. This applies on top of Subject Data Visibility — a sample is shown only when both its subject and the sample itself are marked visible.

Sample_ID	System_Sample_ID	HLA Typing Workshop 2025	Antibody Screening 2025
SAMP_001	IHIW19SAM00000000001	YES	YES
SAMP_002	IHIW19SAM00000000002	YES	YES
SAMP_003	IHIW19SAM00000000003	YES	YES
SAMP_004	IHIW19SAM00000000004	YES	YES

Per-Genotype-sample, per-project visibility matrix

Manage Genotype Sample Visibility

<https://data.ihw19.org/lab/#genotype-sample-visibility>

Lab View

Lab: DOCLAB02

Samples not visible to any project: HLA: 0 / Antibody: 0

Samples visible per project:

- HLA Typing Workshop 2025: HLA 4 (hidden: 0) / Antibody 0 (hidden: 0)
- Antibody Screening 2025: HLA 4 (hidden: 0) / Antibody 0 (hidden: 0)

Lab Members

Project Linkage

Genotype Data

Antibody Testing Data

Sample Metadata

Subject Metadata

Data Visibility

Genotype Sample Visibility

Token Management

Sample_ID / System_Sample_ID: identify each Genotype sample

Control: Individual Genotype samples are visible to each project. Only the sample is shown only when both its subject and the sample itself are visible to the project.

YES (green) = sample visible to project
NO (red) = sample hidden from project

Sample_ID	System_Sample_ID	HLA Typing Workshop 2025	Antibody Screening 2025
SAMP_001	IHIW19SAM00000000001	YES	YES
SAMP_002	IHIW19SAM00000000002	YES	YES
SAMP_003	IHIW19SAM00000000003	YES	YES
SAMP_004	IHIW19SAM00000000004	YES	YES

Manage Genotype Sample Visibility

<https://data.ihw19.org/lab/#genotype-sample-visibility>

IHIW19 [User Guide](#) [Dashboard](#) **[Lab View](#)** [Project View](#)

 Test User ▾

Lab View

Lab: DOCLAB02

Samples not visible to any project: HLA: 0 / Antibody: 0

Samples visible per project:

- HLA Typing Workshop 2025: HLA 4 (hidden: 0) / Antibody 0 (hidden: 0)
- Antibody Screening 2025: HLA 4 (hidden: 0) / Antibody 0 (hidden: 0)


Click "Download Excel"


 [Download Excel](#)

 [Upload Excel](#)

 Lab Members

 Project Linkage

 Genotype Data

 Antibody Testing Data

 Sample Metadata

 Subject Metadata

 Data Visibility

 **Genotype Sample Visibility**

 Token Management

Genotype Sample Visibility

Control which individual Genotype samples are visible to each linked project. This applies on top of Subject Data Visibility — a sample is shown only when both its subject and the sample itself are marked visible.

Sample_ID	System_Sample_ID	HLA Typing Workshop 2025	Antibody Screening 2025
SAMP_001	IHIW19SAM00000000001	YES	YES
SAMP_002	IHIW19SAM00000000002	YES	YES
SAMP_003	IHIW19SAM00000000003	YES	YES
SAMP_004	IHIW19SAM00000000004	YES	YES

Manage Genotype Sample Visibility

<https://data.ihw19.org/lab/#genotype-sample-visibility>

Lab View

Lab: DOCLAB02

Samples not visible to any project: HLA: 0 / Antibody: 0


Samples visible per project:

- HLA Typing Workshop 2025: HLA 4 (hidden: 0) / Antibody 0 (hidden: 0)
- Antibody Screening 2025: HLA 4 (hidden: 0) / Antibody 0 (hidden: 0)

Excel file downloaded successfully

 Lab Members

 Project Linkage

 Genotype Data

 Antibody Testing Data

 Sample Metadata

 Subject Metadata

 Data Visibility

 **Genotype Sample Visibility**

 Token Management

Genotype Sample Visibility

 Download Excel

 Upload Excel

Genotype Sample Visibility Excel downloaded successfully! 

Control which individual Genotype samples are visible to each linked project. This applies on top of Subject Data Visibility — a sample is shown only when both its subject and the sample itself are marked visible.

Sample_ID	System_Sample_ID	HLA Typing Workshop 2025	Antibody Screening 2025
SAMP_001	IHIW19SAM0000000001	YES	YES
SAMP_002	IHIW19SAM0000000002	YES	YES
SAMP_003	IHIW19SAM0000000003	YES	YES
SAMP_004	IHIW19SAM0000000004	YES	YES

Manage HLA Sample Visibility

Open the downloaded Excel file

A	B	C	D	E
Sample_ID	System_Sample_ID	HLA Typing Workshop	Antibody Screening 2025	
SAMP_001	IHIW19SAM0000000001	YES	YES	
SAMP_002	IHIW19SAM0000000002	YES	YES	
SAMP_003	IHIW19SAM0000000003	YES	YES	
SAMP_004	IHIW19SAM0000000004	YES	YES	

“Sample_ID” and “System_Sample_ID” are read-only (gray background)

Manage HLA Sample Visibility

A	B	C	D	E
Sample_ID	System_Sample_ID	HLA Typing Workshop	Antibody Screening 2025	
SAMP_001	IHIW19SAM0000000001	YES	YES	
SAMP_002	IHIW19SAM0000000002	YES	YES	
SAMP_003	IHIW19SAM0000000003	YES	YES	
SAMP_004	IHIW19SAM0000000004	YES	YES	

Project columns: Select YES or NO
YES = HLA sample visible to project leader
NO = HLA sample hidden from project leader

Manage HLA Sample Visibility

A	B	C	D	E
Sample_ID	System_Sample_ID	HLA Typing Workshop	Antibody Screening 2025	
SAMP_001	IHIW19SAM0000000001	NO	NO	
SAMP_002	IHIW19SAM0000000002	YES	YES	
SAMP_003	IHIW19SAM0000000003	YES	NO	
SAMP_004	IHIW19SAM0000000004	YES	YES	

Per-sample, per-project edits:
SAMP_001: hide from BOTH projects
SAMP_003: hide only from Antibody Screening 2025

Only YES and NO values are accepted. Other values will cause an error during upload.

Manage HLA Sample Visibility

A	B	C	D	E
Sample_ID	System_Sample_ID	HLA Typing Workshop	Antibody Screening 2025	
SAMP_001	IHIW19SAM00000000001	NO	NO	
SAMP_002	IHIW19SAM00000000002	YES	YES	
SAMP_003	IHIW19SAM00000000003	YES	NO	
SAMP_004	IHIW19SAM00000000004	YES	YES	

Save the updated Excel file

Manage Genotype Sample Visibility

<https://data.ihw19.org/lab/#genotype-sample-visibility>

Lab View

Lab: DOCLAB02

Samples not visible to any project: HLA: 0 / Antibody: 0

Samples visible per project:

- HLA Typing Workshop 2025: HLA 4 (hidden: 0) / Antibody 0 (hidden: 0)
- Antibody Screening 2025: HLA 4 (hidden: 0) / Antibody 0 (hidden: 0)

Click "Upload Excel"

Download Excel

Upload Excel

Genotype Sample Visibility

Genotype Sample Visibility Excel downloaded successfully!

Control which individual Genotype samples are visible to each linked project. This applies on top of Subject Data Visibility — a sample is shown only when both its subject and the sample itself are marked visible.

Sample_ID	System_Sample_ID	HLA Typing Workshop 2025	Antibody Screening 2025
SAMP_001	IHIW19SAM0000000001	YES	YES
SAMP_002	IHIW19SAM0000000002	YES	YES
SAMP_003	IHIW19SAM0000000003	YES	YES
SAMP_004	IHIW19SAM0000000004	YES	YES

Lab Members

Project Linkage

Genotype Data

Antibody Testing Data

Sample Metadata

Subject Metadata

Data Visibility

Genotype Sample Visibility

Token Management

Manage Genotype Sample Visibility

<https://data.ihw19.org/lab/#genotype-sample-visibility>

Lab View

Lab: DOCLAB02

Samples not visible to any project: HLA: 1 / Antibody: 0


Samples visible per project:

- HLA Typing Workshop 2025: HLA 3 (hidden: 1) / Antibody 0 (hidden: 0)
- Antibody Screening 2025: HLA 2 (hidden: 2) / Antibody 0 (hidden: 0)

Upload completed successfully. The Genotype sample visibility matrix is updated.

 Lab Members

 Project Linkage

 Genotype Data

 Antibody Testing Data

 Sample Metadata

 Subject Metadata

 Data Visibility

 **Genotype Sample Visibility**

 Token Management

Genotype Sample Visibility

 Download Excel

 Upload Excel

Genotype Sample Visibility Excel uploaded successfully! 

Control which individual Genotype samples are visible to each linked project. This applies on top of Subject Data Visibility — a sample is shown only when both its subject and the sample itself are marked visible.

Sample_ID	System_Sample_ID	HLA Typing Workshop 2025	Antibody Screening 2025
SAMP_001	IHIW19SAM0000000001	NO	NO
SAMP_002	IHIW19SAM0000000002	YES	YES
SAMP_003	IHIW19SAM0000000003	YES	NO
SAMP_004	IHIW19SAM0000000004	YES	YES

Manage Genotype Sample Visibility

<https://data.ihiw19.org/lab/#genotype-sample-visibility>

IHIW19 [User Guide](#) [Dashboard](#) **[Lab View](#)** [Project View](#)

 Test User ▾

Lab View

Lab: DOCLAB02


Samples not visible to any project: HLA: 1 / Antibody: 0


Samples visible per project:

- HLA Typing Workshop 2025: HLA 3 (hidden: 1) / Antibody 0 (hidden: 0)
- Antibody Screening 2025: HLA 2 (hidden: 2) / Antibody 0 (hidden: 0)

 Lab Members

 Project Linkage

 Genotype Data

 Antibody Testing Data

 Sample Metadata

 Subject Metadata

 Data Visibility

 **Genotype Sample Visibility**

 Token Management

Genotype Sample Visibility

 Download Excel

 Upload Excel

Visibility settings now apply at the individual Genotype sample level

Control which individual Genotype samples are visible to each linked project on top of Subject Data Visibility — a sample is shown only when both its subject and the sample itself are marked visible.

Sample_ID	System_Sample_ID	HLA Typing Workshop 2025	Antibody Screening 2025
SAMP_001	IHIW19SAM0000000001	NO	NO
SAMP_002	IHIW19SAM0000000002	YES	YES
SAMP_003	IHIW19SAM0000000003	YES	NO
SAMP_004	IHIW19SAM0000000004	YES	YES